DORES BIOLÓGICOS ESTER CONT

Geobacillus Stearothermophilus
(Stearothermophilus) # 7953

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
Geobacillus Stearothermophilus
(BacillusStearothermophilus) # 7953

For: Steam sterilization

Culture: **ESTER CONT** media, 55-59° C. The supplied bacteriological medium will meet requirements for growth promoting ability.

Storage: At controlled room temperature and away from sterilizing agent or other chemicals. Do not desiccate.

Disposal: Incinerate or autoclave at 121°C for not less than 30 minutes.

Purity: No evidence of contaminants using standard plate count techniques.

Lot: 1151UJ114

Expiration: 05/2019

Population: 1,30 . 10⁶ Spores / unit

Assayed Resistance:

Steam Sterilization (121° C ± 1° C)

D-value: 2,14 min. Survival: 8,80 min

Kill: 21,64 min


Z-value: 8,7 °C

D-value and Z-value are reproducible only when exposed and cultured under the exact conditions used to obtain results reported here.

D-value calculated using the Limited- Holcomb-Speaman-Karber method.

Survival and kill values are calculated according to USP e ISO 11138.

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BIOLOGICAL INDICATORS

ESTERCONT Biological Indicators are made up of an external plastic tube with a perforated cap which made it available for the gaseous interchange.

It's contains inside inoculated spore in a filter paper strip and an ampoule with an adequate culture medium and a pH indicator.

This spores are the most adequate for this type of process because they are highly resistant to the sterilization agents and they are not pathogenic for the human been.

To simplify their identification each device have a characteristic colour cap and a label with chemical indicator, which allows to differentiate those that were exposed to the sterilization process of which not were it.

ESTERCONT has been prepared for the O.E and steam sterilization process and allows an exactly biological control with results in a maximum period 48 hours for oxide ethilene and 24 hours for steam.

INSTRUCTION OF USE

1. - EXPOSURE

Place 1 (one) at 4 (four) units of **ESTERCONT** indicators with the material to be sterilized in the most difficult location to sterilize.

Processes the most lead as usual and at the end removed and let cool the indicator by 10 minutes.

Note: The chemical indicator of the label change of yellow to green, when it is exposed to O.E and change of yellow to dark brown when it is exposed to steam. This allows to differentiate those that were exposed to the sterilization process of which not was it.

A colour change doesn't means that the sterilization process was acceptable.

2. - INCUBATION

To activate the media place the indicator in a vertical position, crash the glass ampoule carefully and incubate immediately at 37°C (OE) or 57°C (steam).

3. - INTERPRETATION

Check at regular intervals of time any potential change in colour or turbidity in the culture media (i.e., 12, 18, 24 and 48 hours). The appearance of yellow colour or turbidity indicates bacterial growth. Once finished a period of control of 48 hours for O.E and 24 hours for steam without changes in the culture media can be established a final negative result.

Note: As a positive control and for a lot control, place an indicator no processed and active to incubate.

4. - STORAGE

Store the **ESTERCONT** indicators to room temperature (5°C - 45°C). Don't store near the sterilizing agents or others chemicals.

5. - DISPOSAL

Don't use it after expiration date. Autoclave all positives indicators before eliminating them.